

Congratulations!

You've taken the first step to renewed clarity, confidence and more independence.

Like reading glasses for your ears, our products are designed to enhance your quality of life with a little boost to help you understand what's been missing.

Acclimation

You may not be used to all of the sounds at first, so **taking a little time to get used to them is recommended.**

Start by wearing your Songbird AIR in a quiet area of your house for an hour the first day. Increase your usage to two hours the next, three the following day... and so on.

By the end of the first week, you should be able to comfortably wear your Songbird AIR for six to eight hours a stretch. By the end of week two, you should feel comfortable wearing them from morning to night.

If you have any questions, you can call our customer support team at 1-800-647-5560.

Rest Assured

We're so sure that you'll love your new Songbird Hearing Aids that if you aren't satisfied, just send them back within 45-days of purchase for a full product refund. No questions asked.


Battery Safety

CAUTION – Batteries are dangerous if swallowed. Keep batteries out of reach of children and pets. Never put batteries in your mouth for any reason – they are easily swallowed. In the event a battery has been swallowed, seek immediate medical attention. In the US, in the event a hearing aid or Type 10A battery is swallowed, call collect: The National Button Battery Ingestion Hotline (202) 625-3333 for counsel and treatment.



Songbird AIR

Meet Your New Songbird



Body Hides in your ear and contains both the battery and hearing aid microprocessor.

Battery Door Opens for insertion and removal of the battery. The Battery Door is also the on/off switch. Door open = off. Door closed = on. Your Songbird AIR uses a single size “10A” hearing aid battery.

Tip Sound is emitted from this opening. It is important to keep this sound port clean. You can use our easy to replace screens regularly to keep your Songbird AIR working at its best.

Retention Cord Connects to the body and comfortably rests inside the bowl of your ear to help hold your Songbird AIR in place. If you have larger or smaller ears you may find that using one of the included different-sized retention cords more comfortable.

It is easy to get started

1. Remove your Songbird AIR from the package and attach the appropriate left or right ear retention cord.
2. Insert the battery plus sign (+) up. **When the battery is loaded and the battery door is fully closed, your Songbird is turned on.** To turn it off and conserve the battery, simply open the battery door completely.
3. Inserting your Songbird Air in your ear
 - a). Hold the device with your thumb and index finger while inserting into your ear canal.
 - b). Gently push the retention cord into the back part of the ear.
 - c). Tuck the tip of the retention cord into the top area of your ear and run your finger around the bowl of your ear to “lock” in place.
4. Your Songbird AIR is so advanced that it doesn't require switch or buttons to change volume settings. You can simply **adjust the volume by tapping on your ear two times (2x) in quick succession.** You are trying to generate a “popping” sound, which when heard by your Songbird, will cause it to toggle sequentially to the next setting. **You will hear 1 to 4 beeps as your Songbird AIR changes settings.** One beep = quietest setting. Four beeps = loudest setting. Four loudness settings are preprogrammed for your convenience.

hear here

www.songbirdhearing.com